<b>Romana</b>
<b>Descrierea termometruului</b>
(1) Tasta pornire/oprire
(2) Ecran
(3) Capac baterie
(4) Senzor masurare / varf masurare
(5) Zona de curatat si dezinfectat ( doar sonda )

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**ENGLISH**



**Description of this Thermometer**

- ON/OFF button
- Display
- Battery compartment cover
- Measuring sensor / measuring tip
- Cleaning and disinfecting area (thermometer probe only)


**Important Safety Instructions**

- Follow instructions for use. This document provides important product operation and safety information regarding this device.
- Please read this document thoroughly before using the device and keep for future reference.
- This device is only to be used for measuring human body temperature through oral, rectal or axillary. Do not attempt to take temperatures at other sites, such as in the ear, as it may result in false readings and may lead to injury.
- Do not use this device if you think it is damaged or notice anything unusual.
- We recommend cleaning this device according to the cleaning instructions before first use for personal hygiene.
- The minimum measurement time until the beep is heard must be maintained without exception! Consider that different measurement locations may require continued measuring even after the beep, see section «Measuring methods / Normal body temperatures».
- For safety reasons (risk of rectal perforation) rectal measurement in children younger than 3 years must be performed only by trained health care personnel (professional user). Use another measurement method instead. For rectal fever measurement in small children younger than 3 years, thermometers with a flexible tip are available.
- Do not attempt rectal measurements on persons with rectal disorders.
- Doing so may aggravate or worsen the disorder.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect the device from impact and dropping!
- Avoid ambient temperatures above 60 °C. NEVER boil this device!
- Use only the commercial disinfectants listed in the section «Cleaning and Disinfecting» to clean the device to avoid damage to the device.
- We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact your Riester-Service to arrange the test.

<span><span>⚠</span></span>	<b>WARNING:</b> The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only.
<span><span>🗑</span></span>	Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.
<span><span>👤</span></span>	Read the instructions carefully before using this device.
<span><span>👤</span></span>	Type BF applied part

- (1) ON/OFF button**  
**(2) Display**  
**(3) Battery compartment cover**  
**(4) Measuring sensor / measuring tip**  
**(5) Cleaning and disinfecting area (thermometer probe only)**

**CE 0044**

🇪🇨 REP **Microlife AG**  
Espanenstrasse 139  
9443 Widnau  
Switzerland

Distributed by / Vertrieb durch:  
**Rudolf Riester GmbH**  
Bruckstr. 31  
72471 Jungingen | Germany  
info@riester.de  
www.riester.de

**Turning on the Thermometer**
To turn on the thermometer, press the ON/OFF button (1); a short beep signals «thermometer ON». A display test is performed. All segments should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «°C» appear at the display field (2). The thermometer is now ready for use.

**Function Test**

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

**Using the Thermometer**

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «°C» symbol flashes. If the beep is heard 10 times and the «°C» is no longer flashing, this means that the measured increase in temperature is less than 0.1 °C in 16 seconds. To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button (1). Otherwise the thermometer will automatically turn off after about 10 minutes.

**Storage of Measured Values**

If the ON/OFF button (1) is pressed for more than 3 seconds when turning on the thermometer, the automatically stored maximum temperature during the last measurement will be displayed. At the same time, a «M» for memory will appear on the display. About 2 seconds after the button is released, the temperature value disappears and the thermometer is ready for measurement.

**Measuring methods / Normal body temperature**

**În the armpit (axillary) / 34,7 - 37,3 °C**
Wipe the underarm with a dry towel. Place the measuring sensor (4) under the arm into the center of the armpit so the tip is touching the skin and position the patient's arm next to the patient's body. This ensures that the room air does not affect the reading. Because the axillary takes more time to reach ist stable temperature wait **at least 5 minutes**, regardless of the beep sound.

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Do not eat or drink anything hot or cold 10 minutes before the measurement. The mouth should remain closed up to 2 minutes before starting a reading. Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor (4) must be in good contact with the tissue. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. If this is not possible due to blocked airways, another method for measuring should be used.

**Approx. measuring time: 1 minute!**

**În the anus (rectal) / 36.6 - 38.0 °C**
**Attention:** For the prevention of rectal perforation in children (younger than 3 years), we recommend using another measuring method, or the use of a thermometer with flexible tip. Carefully insert the measuring sensor (4) of the thermometer 2 to 3 cm into the anal aperture. The use of a probe cover and the use of a lubricant is recommended. If you are unsure of this measurement method, you should consult a professional for guidance/training.
**Approx. measuring time: 1 minute!**

**Cleaning and Disinfecting**

For disinfection in home use environment, use a 70% Isopropyl alcohol swab, or a cotton tissue moistened with 70% Isopropyl alcohol to wipe surface pollutants off the thermometer probe (note: consider the application and safety instruction of the disinfectant manufacturer). Always start wiping from the end of the thermometer probe (approx. at the middle of the thermometer) towards the thermometer tip. Afterwards the entire thermometer probe (see number (5) in the drawing) should be immersed in 70% Isopropyl alcohol for at least 5 minutes (max. 24 hours).After immersion, let the disinfectant dry off for 1 minute before next use. Avoid immersing or wiping the display to protect it from fading.

The thermometer is not intended for professional use.

**Battery Replacement**

When the «▼» symbol (upside-down triangle) appears in the display, the battery is flat and needs replacing. To replace the battery remove the battery compartment cover (3) from the thermometer. Insert the new battery with the + at the top. Make sure you have a battery of the same type to hand. Batteries can be purchased at any electrical store.

<b>Technical Specifications</b>	
<b>Type:</b>	Maximum thermometer
<b>Measurement range:</b>	32,0 <span> </span> °C to 43,9 <span> </span> °C
	Temp. < 32,0 <span> </span> °C: display «L» for low (too low)
	Temp. > 43,9 <span> </span> °C: display «H» for high (too high)
	± 0.1 <span> </span> °C; 34 <span> </span> °C - 42 <span> </span> °C
	± 0.2 <span> </span> °C; 32,0 - 33,9 <span> </span> °C and 42.1 - 43,9 <span> </span> °C
	10 - 40 <span> </span> °C; 15 - 95 <span> </span> % relative maximum humidity
<b>Operating conditions:</b>	-25 - +60 <span> </span> °C; 15 - 95 <span> </span> % relative maximum humidity
<b>Storage conditions:</b>	LR41 [1.5 V] / SR41 [1.55 V]
<b>Battery:</b>	approx. 4.500 measurements (using a new battery)
<b>Battery lifetime:</b>	IP22
<b>IP Class:</b>	
<b>Reference to standards:</b>	EN 12470-3, clinical thermometers; ASTM E1112; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11
<b>Expected service life:</b>	5 years or 10.000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

**Warranty**

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of **2 years from the date of purchase** on all defects, which can verifiably be shown to be due to material or manufacturing faults. A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts. A warranty claim can only be granted if this Warranty Card has been completed and stamped by the dealer and is enclosed with the product. Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the RIESTER product with the completed Warranty Card to the following address:

**Rudolf Riester GmbH**  
**Dept. Repairs RR**  
**Bruckstr. 31**  
**72471 Jungingen**  
**Germany**

**Serial number or batch number / Date / Stamp and signature of the specialist dealer**

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