

Traducere eticheta / instructiuni utilizare

Ultraseptic

ENGLEZA: ULTRASOUND PROBE CLEANING SOLUTION

INTENDED USE: Enables maximum visibility due to the cleaning gel residues on the probe and extends the service life of the probe.

USAGE: It is sprayed onto the probe and wiped with a clean cloth or tissue paper before and after the ultrasound procedure.

INGREDIENTS: Deionized water, Nonylphenol Ethoxylate, Benzalkonium Chloride, Methylchloroisothiazolineone, Methylisothiazolineone

ROMANA : SOLUTIE DE CURATARE SONDE DE ULTRASUNET

Scop utilizare recomandat: Permite o vizibilitate maximă datorită reziduurilor de gel de curățare de pe sondă și prelungeste durata de viață a sondei.

Mod utilizare: Se pulverizează pe sondă și se șterge cu o cârpă curată sau cu un șervețel de hârtie înainte și după procedura cu ultrasunete.

INGREDIENTE: Apă deionizată, Etoxilat de nonilfenol, Clorură de benzalconiu, Metilcloroizotiazolinonă, Metilizotiazolinonă.

